

GENADYNE BIOTECHNOLOGIES

PRODUCT CATALOG

Rev A.1.1 October 2025

GENADYNE

The power to *Heal*

TABLE OF CONTENT

<hr/>		<u>Multi-patient use NPWT.</u>		<hr/>	
	<u>DUO</u>				
Therapy unit			Page		3.
Canisters			Page		3.
Dressing kits			Page		4-7.
Individual dressings			Page		8-9.
			Page		10.
Accessories					
<hr/>		<u>Multi-patient use NPWT.</u>		<hr/>	
	<u>XLR8+</u>				
Therapy unit			Page		11.
Canisters			Page		12.
Dressing kits			Page		13-16.
Individual dressings			Page		17-18.
Accessories			Page		19.
<hr/>		<u>Multi-patient use NPWT.</u>		<hr/>	
	<u>UNO+</u>				
Therapy unit			Page		20.
Canisters			Page		20.
Dressing kits			Page		21-24.
Individual dressings			Page		24-25.
Accessories			Page		26.
<hr/>		<u>Single-patient use NPWT.</u>		<hr/>	
	<u>UNO</u>				
Therapy unit			Page		27.
Canisters			Page		28.
Dressing kits			Page		29-30.
<hr/>		<u>Advanced wound care.</u>		<hr/>	
Gentilfoam- Gentiltac			Page		31-32.
Nanogen			Page		32.
Silgen			Page		33.

DUO



Duo NPWT-id system.
Ref. Code: A7-S0001.



- 3 Intensity Setting Options.
- Instill, Continuous or Variable Pressure.
- Lightweight:2.60kg.
- Small size: 178.5 x 243 x 168 mm.
- Whisper Quiet Operation.
- Maximum recharge time of 2.5 hours
- Battery life up to 10 hours
- Pressure range:
 1. Continuous 25 mmHg – 200mmHg.
 2. Variable: 0 mmHg - 200mm Hg.
- Supports Multiple Languages.
- 8 Advanced Alert Notifications: Leakage, blockage, canister Full, low battery, critical battery, solution bag empty, instill line blockage, instill door open.
- Clear and high visibility display.
- Simple, lockable controls.

The Genadyne DUO Negative Pressure Wound Therapy (NPWT) system combines advanced wound care technology with instillation therapy, allowing liquid solutions to be delivered directly to the wound bed and set to dwell for optimal cleansing. With continuous and variable pressure modes, DUO supports precision-driven healing tailored to patient needs.

CANISTERS



Duo canister 600cc.
Ref. Code: A7-S00D6.



Duo canister 800cc.
Ref. Code: A7-S00D8.



Duo canister 1100cc.
Ref. Code: A7-S00D11.

DRESSING KITS

DUO



DUO INSTILL NPWT dressing kit.

Comes with instill foam, 1 instill tubing and 1 port pad (Dome diameter 5cm and disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Foam	Size cm	Kit / Case	Film /Kit	Picture
DIF-SPSMF1	Spiro instill	Small (105 x 2 x 1.5)	10 / Case	2	
DIF-SPMDF1	2 Spiro instill	Small (105 x 2 x 1.5)	10 / Case	3	
DIF-SPLGF1	Large foam instill	Large (23 x 15 x 1.5)	10 / Case	5	
DIF-HFLGF1	Large foam instill with holes	2 Large thin (26 x 15 x 1) 1 Large thin with holes (26 x 15 x 1)	10 / Case	5	

DRESSING KITS

DUO



DUO NPWT dressing kit

Comes with 1 green foam and 1 port pad (Dome diameter 5cm and disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Size cm	Kit / Case	Film /Kit	Picture
XF-DSMF1-D	Small (10 x 7.5 x 3.3)	10 / Case	1	
XF-DMDF1-D	Medium (18 x 12.5 x 3.3)	10 / Case	2	
XF-DLGF1-D	Large (26 x 15 x 3.3)	10 / Case	2	
XF-DXLF1-T-D	Thin extra large (60 x 30 x 1.5)	10 / Case	5	

DUO NPWT dressing Kit.

Comes with 1 green foam and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Size (cm)	Kit / Case	Film / Kit	Picture
XF-DSMF1-T-D	Thin small (6 x 6 x 2.5)	10 / Case	1	
XF-DMDF1-T-D	Thin medium (15 x 10 x 2)	10 / Case	2	
XF-DLGF1-T-D	Thin large (23 x 15 x 2.5)	10 / Case	2	

DRESSING KITS

DUO



Green spiro foam dressing kit.

Comes with foam and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Size (cm)	Foams	Kit / Case	Film / Kit	Picture
XF-SPSK1-D	Spiro foam small (105 x 2 x 1.5)	1	10 / Case	1	
XF-SPMK1-D	Spiro foam medium 2 x (105 x 2 x 1.5)	2	10 / Case	2	

DUO PVA dressing Kit.

Comes with 1 PVA foam and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Size (cm)	Kit / Case	Film / Kit	Picture
PVA-SFOAM1-D	Small (10 x 7.5 x 1)	10 / Case	1	
PVA-FOAM1-D	Medium (15 x 10 x 1)	10 / Case	2	
PVA-LFOAM1-D	Large (20 x 15 x 1)	10 / Case	2	

DUO bridge dressing kit.

Comes with 1 bridge, 1 foam and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Size (cm)	Kit / Case	Film / Kit	Bridge
XF-BKSF1-D	Small (15 x 10 x 2)	10 / Case	1	 Length 30 cm
XF-BKMF1-D	Medium (15 x 10 x 2)	10 / Case	1	 Length 55 cm


DRESSING KITS

DUO



DUO Abdominal dressing kit.

Comes with 2 foams, 1 Y connector and 2 port pads (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm) and a visceral protective layer (short axis 61 x long axis 79).

Ref. Code	Size (cm)	Kit / Case	Film / Kit	Bridge (Length 55cm)
A4-ABDO1-K-D	2 Oval shaped foam (39.1 x 26.1 x 1.7)	10 / Case	4	

DUO silver dressing Kit.

Comes with 1 foam and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Size (cm)	Kit / Case	Film / Kit	Picture
AG-SFOAM2-D	Small (7.5 x 10 x 3.3)	10 / Case	1	
AG-MFOAM2-D	Medium (12.5 x 18 x 3.3)	10 / Case	2	
AG-LFOAM2-D	Thin large (15 x 26 x 3.3)	10 / Case	2	
XF-SPSK1-AG-D	Spiro foam 1 x (105 x 2 x 1.5)	10 / Case	1	
XF-SPMK1-AG-D	Spiro foam 1 x (105 x 2 x 1.5)	10 / Case	2	

INDIVIDUAL DRESSINGS

DUO



White PVA foam

Ref. Code	Description and size (cm)	Qty / Case	Picture
XF-PVASFOAM1	PVA small foam (10 x 7.5 x 1)	10 / Case	
XF-PVAFOAM1	PVA medium foam (15 x 10 x 1)	10 / Case	
XF-PVALFOAM1	PVA large foam (20 x 15 x 1)	10 / Case	

Green spiro foam.

Ref. Code	Size (cm)	Foams	Kit / Case	Picture
XF-SPSF1	Spiro foam small (105 x 2 x 1.5)	1	10 / Case	
XF-SPMF1	Spiro foam medium 2 x (105 x 2 x 1.5)	2	10 / Case	

Green thick round foam.

Ref. Code	Size (cm)	Kit / Case	Picture
XF-SFOAM2	Small (7.5 x 10 x 3.3)	10 / Case	
XF-MFOAM2	Medium (12.5 x 18 x 3.3)	10 / Case	
XF-LFOAM2	Large (15 x 26 x 3.3)	10 / Case	

INDIVIDUAL DRESSINGS

DUO



Green thin square foam.

Ref. Code	Size (cm)	Unit / Case	Picture
XF-SFOAM1	Small (6 x 6 x 2.5)	10 / Case	
XF-MFOAM1	Medium (15 x 10 x 2)	10 / Case	
XF-LFOAM1	Large (15 x 23 x 2.5)	10 / Case	
XF-XLFOAM1	(30 x 60 x 1.5)	10 / Case	

XLR8+ Bridge dressing.

Ref. Code	Size (cm)	Unit/ Case	Picture
XF-BDSF1	Length 30	10 / Case.	
XF-BDMF1	Length 55	10 / Case.	

ACCESSORIES

DUO



DUO Accessories

Ref. code	Description	Qty	Picture
DP-1000	DUO port pad double lumen	10 / Case	
DP-2000	DUO port pad triple lumen	10 / Case	
DP-3000	DUO instill tubing	10 / Case	
A7-S00Y2	DUO NPWT Y connector	10 / Case	
A4-S00F5	Transparent film dressing (26 x 30)	100 / Case.	

XLR8+

XLR8+ NPWT system.
Ref. code: A4-S0003



The Genadyne XLR8+ takes NPWT to the next level with enhanced power and performance. Engineered to handle the largest wounds with ease, it offers continuous and variable therapy along with intensity level settings for extended treatment durations.



Specifications	
Therapy unit	XLR8+ NPWT
Weight	< 1.65 lbs
Suction capacity	~ 5.2 to 6.4 Liters/Minute
Continuous & variable therapy	✓
Whisper quiet operation	✓
Supports multiple languages	✓
5 Advanced alerts notifications	✓
Rechargeable lithium battery	✓
Lockable controls	✓
Global regulatory approved	✓
Intensity controls	✓
Upgradeable for wireless capability	✓

CANISTERS

XLR8+








Multi-capacity canister system

Allows the clinician to choose the appropriate size canister for the amount of wound drainage.

Secure luer lock attachment to port pad.

Mix and match capability – use any canister with any dressing kit.

All sizes compatible with carrying bag.

Ref. code	Volume	Picture
A4-S00D2	200cc	
A4-S00D4	400cc	
A4-S00D6	600cc	
A4-S00D8	800cc	
A4-S00D11	1100cc	

DRESSING KITS

XLR8+



Green foam dressing Kit.

Comes with 1 foam and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Size (cm)	Kit / Case	Film / Kit	Picture
XF-DSMF1	Small (7.5 x 10 x 3.3)	10 / Case	1	
XF-DMDF1	Medium (12.5 x 18 x 3.3)	10 / Case	2	
XF-DLGF1	Large (15 x 26 x 3.3)	10 / Case	2	
XF-DXLF1-T	Thin extra large (30 x 60 x 1.5)	10 / Case	5	
XF-DSMF1-T	Thin small (6 x 6 x 2.5)	10 / Case	1	
XF-DMDF1-T	Thin medium (10 x 15 x 2)	10 / Case	2	
XF-DLGF1-T	Thin large (15 x 23 x 2.5)	10 / Case	2	

Green spiro foam dressing kit.

Comes with foam and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Size (cm)	Foams	Kit / Case	Film / Kit	Picture
XF-SPSK1	Spiro foam small (105 x 2 x 1.5)	1	10 / Case	1	
XF-SPMK1	Spiro foam medium 2 x (105 x 2 x 1.5)	2	10 / Case	2	

DRESSING KITS

XLR8+



PVA foam dressing Kit.

Comes with 1 foam and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Size (cm)	Kit / Case	Film / Kit	Picture
PVA-SFOAM1	Small (10 x 7.5 x 1)	10 / Case	1	
PVA-FOAM1	Medium (15 x 10 x 1)	10 / Case	2	
PVA-LFOAM1	Thin large (20 x 15 x 1)	10 / Case	2	

Bridge dressing kit.

Comes with 1 bridge, 1 foam and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Size (cm)	Kit / Case	Film / Kit	Bridge
XF-BKSF1	Small (15 x 10 x 2)	10 / Case	1	 Length 30 cm.
XF-BKMF1	Medium (15 x 10 x 2)	10 / Case	1	 Length 55 cm.

DRESSING KITS

XLR8+



XLR8+ Silver foam dressing Kit.

Comes with foam and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Size (cm)	Kit / Case	Film / Kit	Picture
AG-SFOAM2	Small foam (7.5 x 10 x 3.3)	10 / Case	1	
AG-MFOAM2	Medium foam (12.5 x 18 x 3.3)	10 / Case	2	
AG-LFOAM2	Large foam (15 x 26 x 3.3)	10 / Case	2	
XF-SPSK1-AG	Spiro foam (104 x 2 x 1.5)	10 / Case	1	
XF-SPMK1-AG	Spiro foam 2 x (104 x 2 x 1.5)	10 / Case	2	

XLR8+ Abdominal dressing kit.

Comes with 2 foams, 1 Y connector and 2 port pads (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm) and a visceral protective layer (short axis 61 x long axis 79).

Ref. Code	Size (cm)	Kit / Case	Film / Kit	Picture
A4-ABD01-K	2 Oval shaped foam (39.1 x 26.1 x 1.7)	5 / Case	4	

DRESSING KITS

XLR8+



XLR8+ Gauze foam dressing Kit.

Comes with 1 gauze and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm), 1 protective barrier and oil emulsion.

Ref. Code	Size (cm)	Kit / Case	Oil emulsion Qty & size	Film / Kit	Picture
A4-D00SP	Small gauze (15 x 17.1)	10 / Case	1 (7.6 x 7.6)	1	
A4-D00MP	Medium gauze (15 x 17.1)	10 / Case	2 (7.6 x 7.6)	2	
A4-D00LP	Large gauze (10.16 x 365)	10 / Case	2 (7.6 x 23)	2	

Hybrid Dressing

comes with 1 hybrid dressing, 1 port pad and 5 film strips.

Ref. Code	Size (cm)	Kit / Case	Picture
UP-D1020	10 x 20	5 / Case	
UP-D1030	10 x 30	5 / Case	
UP-D1040	10 x 40	5 / Case	
UP-D1520	15 x 20	5 / Case	
UP-D1530	15 x 30	5 / Case	
UP-D1050	10 x 50	5 / Case	

INDIVIDUAL DRESSINGS

XLR8+



White PVA foam

Ref. Code	Description and size (cm)	Qty / Case	Picture
XF-PVASFOAM1	PVA small foam (10 x 7.5 x 1)	10 / Case	
XF-PVAFOAM1	PVA medium foam (15 x 10 x 1)	10 / Case	
XF-PVALFOAM1	PVA large foam (20 x 15 x 1)	10 / Case	

Green spiro foam.

Ref. Code	Size (cm)	Foams	Kit / Case	Picture
XF-SPSF1	Spiro foam small (105 x 2 x 1.5)	1	10 / Case	
XF-SPMF1	Spiro foam medium 2 x (105 x 2 x 1.5)	2	10 / Case	

Green thick round foam.

Ref. Code	Size (cm)	Kit / Case	Picture
XF-SFOAM2	Small (7.5 x 10 x 3.3)	10 / Case	
XF-MFOAM2	Medium (12.5 x 18 x 3.3)	10 / Case	
XF-LFOAM2	Large (15 x 26 x 3.3)	10 / Case	

INDIVIDUAL DRESSINGS

XLR8+



Green thin square foam.

Ref. Code	Size (cm)	Unit / Case	Picture
XF-SFOAM1	Small (6 x 6 x 2.5)	10 / Case	
XF-MFOAM1	Medium (15 x 10 x 2)	10 / Case	
XF-LFOAM1	Large (15 x 23 x 2.5)	10 / Case	
XF-XLFOAM1	(30 x 60 x 1.5)	10 / Case	

XLR8+ Bridge dressing.

Ref. Code	Size (cm)	Unit/ Case	Picture
XF-BDSF1	Length 30	10 / Case	
XF-BDMF1	Length 55	10 / Case	

ACCESSORIES

XLR8+



Ref. Code	Description and size (cm)	Qty / Case	Picture
A4-S00F5	Transparent film dressing (26 x 30)	100 / Case	A rectangular, transparent adhesive film with a grid pattern of small red markings.
A4-S00A9	XLR8+ bed hanger	1	A white, wire-mesh basket with a handle on top, designed for hanging a device.
A4-S00A7	XLR8+ Carrying bag	1	A black, soft carrying bag with a shoulder strap and a handle. It contains the device and some clear tubing.
A4-S00A10	XLR8+ cordura bag	1	A black, soft carrying bag with a shoulder strap and a handle, similar to the previous one but without the device inside.
X-A0Y1-P	Y connector	10 / Case	A small, clear plastic Y-shaped connector with two ports on one side and one on the other.
XP-1013	XLR8+ flat port pad	10 / Case	A small, circular, flat port pad with a clear tube attached to one side.

UNO+



UNO Plus NPWT system.
Ref. Code: UP-S01.



The Genadyne UNO+ is a compact and versatile Negative Pressure Wound Therapy (NPWT) device designed for simplicity and efficiency. It delivers continuous and variable therapy while offering a lightweight, discreet design that enhances patient mobility and comfort

Discreet, comfortable and lightweight.

- 220 grams or 0.49 lbs.
- Small and portable size of 111.13 mm x 76.2 mm x 41.28 mm or 4 3/8 " x 3 " x 1 5/8".
- Proprietary PC Interface for downloading alert log and other diagnostic data.
- Up to 40 hours of usage per charge, max recharge time of 3 hours.
- Two therapy settings:
 1. Continuous at 40mmHg, 60mmHg, 80mmHg, 100mmHg and 125mmHg
 2. Variable at 80mmHg / 30mmHg, 100mmHg / 30mmHg or 125mmHg / 30mmHg
- 4 Alerts: Leakage ,Canister Full, Blockage, Low Battery.

CANISTERS

UNO Plus 100cc canister.
Ref. Code: UP-S00D1.



UNO Plus 200cc canister.
Ref. Code: UP-S00D2.



UNO Plus 300cc canister.
Ref. Code: UP-S00D3.



ACCESSORIES

UNO Plus carrying bag.
Ref. Code: UP-CBAG.



DRESSING KITS



Green foam dressing Kit.

Comes with 1 foam and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Size (cm)	Kit / Case	Film / Kit	Picture
XF-DSMF1	Small (7.5 x 10 x 3.3)	10 / Case	1	
XF-DMDF1	Medium (12.5 x 18 x 3.3)	10 / Case	2	
XF-DSMF1-T	Thin small (6 x 6 x 2.5)	10 / Case	1	
XF-DMDF1-T	Thin medium (10 x 15 x 2)	10 / Case	2	

Green spiro foam dressing kit.

Comes with foam and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Size (cm)	Foams	Kit / Case	Film / Kit	Picture
XF-SPSK1	Spiro foam small (105 x 2 x 1.5)	1	10 / Case	1	
XF-SPMK1	Spiro foam medium 2 x (105 x 2 x 1.5)	2	10 / Case	2	

DRESSING KITS

UNO+



UNO+ Silver foam dressing Kit.

Comes with foam and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Size (cm)	Kit / Case	Film / Kit	Picture
AG-SFOAM2	Small foam (7.5 x 10 x3.3)	10 / Case	1	
AG-MFOAM2	Medium foam (12.5 x 18 x 3.3)	10 / Case	2	
XF-SPMK1-AG	Spiro foam 2 x (104 x 2 x 1.5)	10 / Case	2	

Hybrid Dressing

comes with 1 hybrid dressing, 1 port pad and 5 film strips.

Ref. Code	Size (cm)	Kit / Case	Picture
UP-D1020	10 x 20	5 / Case	
UP-D1030	10 x 30	5 / Case	
UP-D1040	10 x 40	5 / Case	
UP-D1520	15 x 20	5 / Case	
UP-D1530	15 x 30	5 / Case	
UP-D1050	10 x 50	5 / Case	

DRESSING KITS



PVA foam dressing Kit.

Comes with 1 foam and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Size (cm)	Kit / Case	Film / Kit	Picture
PVA-SFOAM1	Small (10 x 7.5 x 1)	10 / Case	1	
PVA-FOAM1	Medium (15 x 10 x 1)	10 / Case	2	

Bridge dressing kit.

Comes with 1 bridge, 1 foam and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm).

Ref. Code	Size (cm)	Kit / Case	Film / Kit	Bridge
XF-BKSF1	Medium (15 x 10 x 2)	10 / Case	1	 Length 30 cm.
XF-BKMF1	Medium (15 x 10 x 2)	10 / Case	1	 Length 55 cm.

DRESSING KITS

UNO+



UNO+ Gauze foam dressing Kit.

Comes with 1 gauze and 1 port pad (Dome diameter 5 cm, Disc film 10 x 13 cm, Tube length 80 cm), 1 protective barrier and oil emulsion.

Ref. Code	Size (cm)	Kit / Case	Oil emulsion Qty & size	Film / Kit	Picture
A4-D00SP	Small gauze (15 x 17.1)	10 / Case	1 (7.6 x 7.6)	1	
A4-D00MP	Medium gauze (15 x 17.1)	10 / Case	2 (7.6 x 7.6)	2	

INDIVIDUAL DRESSINGS

Green spiro foam.

Ref. Code	Size (cm)	Foams	Kit / Case	Picture
XF-SPSF1	Spiro foam small (105 x 2 x 1.5)	1	10 / Case	
XF-SPMF1	Spiro foam medium 2 x (105 x 2 x 1.5)	2	10 / Case	

Green thick round foam.

Ref. Code	Size (cm)	Kit / Case	Picture
XF-SFOAM2	Small (7.5 x 10 x 3.3)	10 / Case	
XF-MFOAM2	Medium (12.5 x 18 x 3.3)	10 / Case	

INDIVIDUAL DRESSINGS



White PVA foam

Ref. Code	Description and size (cm)	Qty / Case	Picture
XF-PVASFOAM1	PVA small foam (10 x 7.5 x 1)	10 / Case	
XF-PVAFOAM1	PVA medium foam (15 x 10 x 1)	10 / Case	

Green thin square foam.

Ref. Code	Size (cm)	Unit / Case	Picture
XF-SFOAM1	Small (6 x 6 x 2.5)	10 / Case	
XF-MFOAM1	Medium (15 x 10 x 2)	10 / Case	

UNO+ Bridge dressing.

Ref. Code	Size (cm)	unit/ Case	Picture
XF-BDSF1	Length 30	10 / Case.	
XF-BDMF1	Length 55	10 / Case.	

ACCESSORIES

UNO+



Ref. Code	Description and size (cm)	Qty / Case	Picture
X-A0Y1-P	Y connector	10 / Case	
XP-1013	XLR8 flat port pad	10 / Case	

UNO

UNO single-patient NPWT system.

Ref. Code: U-S00R4/ U-S00R5/
U-S00R6



Genadyne's UNO single-patient NPWT system offers a rechargeable Therapy unit with continuous and variable therapy. UNO manages and distributes excess exudate, through the dressing and into the canister, ensuring most levels of exudate are managed with ease.

Discreet, comfortable and lightweight.




- 7, 15 and 30 day operating options available.
- Discreet, comfortable and lightweight: 230 grams or 0.5 lbs
- Usage of ~ 48-72 hours per full charge.*
- Simple, easy, clear and intuitive interface.
- Two therapy setting options:
 1. Continuous at 80mmHg or 125mmHg.
 2. Variable at 80mmHg / 30mmHg or 125mmHg / 30mmHg.
- 3 Alerts: Leakage, canister full, low battery.

Ref. code	Description	Picture
U-S00R4	7 days UNO NPWT Therapy unit	
U-S00R5	15 days UNO NPWT Therapy unit	
U-S00R6	30 days UNO NPWT Therapy unit	

PRICE

UNO

CANISTER



UNO canister 70ml set of 2.
Ref. Code: U-C070




ACCESSORIES

UNO carry pouch.
Ref. Code: U-000A2-R



DRESSING



Polyurethane film layer.

Allows small amounts of air to pass through, needed for exudate to be pulled through into the canister. Transpires exudate away from the wound bed.

Super absorbant layer.

For additional absorption and aids the movement of exudate into the port.

Foam layer.

For quick absorption and maintaining a moist environment at the wound site.

Silicone wound contact layer.

Soft and gentle application and removal for enhanced patient comfort. Can be repositioned to guarantee a good seal on the dressing.

THERAPY KITS

UNO



UNO 7 therapy kit

(Comes with 1 UNO 7 day Therapy unit, charger and carry case.)



Ref. code	Size	Kit trays	Foam kit/Tray	Canister 70mL/Case
U-S1010-4	10x10	2 / Case	1 Set	1
U-S1020-4	10x20	2 / Case	1 Set	1
U-S1030-4	10x30	2 / Case	1 Set	1
U-S1040-4	10x40	2 / Case	1 Set	1
U-S1515-4	15x15	2 / Case	1 Set	1
U-S1520-4	15x20	2 / Case	1 Set	1
U-S1530-4	15x30	2 / Case	1 Set	1
U-S2020-4	20x20	2 / Case	1 Set	1
U-S2025-4	20x25	2 / Case	1 Set	1
U-S2525-4	25x25	2 / Case	1 Set	1

UNO 15 therapy kit

(Comes with 1 UNO 15 day Therapy unit, charger and carry case.)



Ref. code	Size	Kit trays	Foam kit/Tray	Canister 70mL/Case
U-S1010-5	10x10	5 / Case	1 Set	2
U-S1020-5	10x20	5 / Case	1 Set	2
U-S1030-5	10x30	5 / Case	1 Set	2
U-S1040-5	10x40	5 / Case	1 Set	2
U-S1515-5	15x15	5 / Case	1 Set	2
U-S1520-5	15x20	5 / Case	1 Set	2
U-S1530-5	15x30	5 / Case	1 Set	2
U-S2020-5	20x20	5 / Case	1 Set	2
U-S2025-5	20x25	5 / Case	1 Set	2
U-S2525-5	25x25	5 / Case	1 Set	2

THERAPY KITS

UNO



UNO 30 therapy kit **i**

(Comes with 1 UNO 30 day Therapy unit, charger and carry case.)

Ref. code	Size	Kit trays	Foam kit/Tray	Canister 70mL/Case
U-S1010-6	10x10	5 / Case	1 Set	2
U-S1020-6	10x20	5 / Case	1 Set	2
U-S1030-6	10x30	5 / Case	1 Set	2
U-S1040-6	10x40	5 / Case	1 Set	2
U-S1515-6	15x15	5 / Case	1 Set	2
U-S1520-6	15x20	5 / Case	1 Set	2
U-S1530-6	15x30	5 / Case	1 Set	2
U-S2020-6	20x20	5 / Case	1 Set	2
U-S2025-6	20x25	5 / Case	1 Set	2
U-S2525-6	25x25	5 / Case	1 Set	2


UNO dressing kits **i**

Ref. code	Size	Kit trays	Foam kit/Tray	Film strips
UNO-0404	10x10	5 / Case	1	6
UNO-0408	10x20	5 / Case	1	6
UNO-0412	10x30	5 / Case	1	6
UNO-0416	10x40	5 / Case	1	6
UNO-0606	15x15	5 / Case	1	6
UNO-0608	15x20	5 / Case	1	6
UNO-0612	15x30	5 / Case	1	6
UNO-0808	20x20	5 / Case	1	6
UNO-0810	20x25	5 / Case	1	6
UNO-1010	25x25	5 / Case	1	6




AWC PRODUCTS

Gentil foam with border Soft Silicone foam dressing.

Ref. Code	Size (cm)	Qty / Carton	Picture
430303	Border (7.5 x 7.5)	10 / Carton	
430404	Border (10 x 10)	10 / Carton	
430505	Border (12.5 x 12.5)	10 / Carton	
430606	Border (15 x 15)	10 / Carton	
430608	Border (15 x 20)	10 / Carton	


Gentil foam without border Soft Silicone foam dressing.

Ref. Code	Size (cm)	Qty / Carton	Picture
420202	5 x 5	10 / Carton	
420303	7.5 x 7.5	10 / Carton	
420404	10 x 10	10 / Carton	
420606	15 x 15	10 / Carton	
420608	15 x 20	5 / Carton	

AWC PRODUCTS

GentilTac

Soft silicone interface layer.



Ref. Code	Description and size (cm)	Qty / Carton	Picture
440203	5 x 7.5	10 / Carton	
440304	7.5 x 10	10 / Carton	
440407	10 x 18	10 / Carton	
440610	17 x 25	10 / Carton	
440812	20 x 30	5 / Carton	
440911	24 x 27.5	5 / Carton	

Nanogen Aktiv membrane

Bio-cellulose nano-structured matrix.

Ref. Code	Description and size (cm) or Volume.	Qty / Carton	Picture
610101	Membrane (3 x 3)	10 / Carton	
610202	Membrane (5 x 5)	10 / Carton	
610303	Membrane (8 x 8)	5 / Carton	
610404	Membrane (10 x 10)	5 / Carton	
610408	Membrane (10 x 20)	5 / Carton	
620015	Aktigel 15mL	5 / Carton	
620090	Aktigel 90 mL	5 / Carton	

AWC PRODUCTS



Silgen Silver spray

Ref. Code	Description Volume.	Qty / Carton	Picture
710005	Spray bottle 5 mL.	10 / Carton	
710010	Spray bottle 10 mL.	5 / Carton	
710025	Spray bottle 25 mL.	5 / Carton	
710100	Spray bottle 100 mL.	1 / Carton	
S-710005	Spray bottle 5 mL in sterile bag	5 / Carton	
S-710010	Spray bottle 10 mL in sterile bag	5 / Carton	



Genadyne Biotechnologies, Inc.

16 Midland Ave., Hicksville, NY 11801

www.genadyne.com

If you have questions about our product or need support, or interested in becoming a distributor, please email info@genadyne.com or call +1.516.487.8787 / +1.800.208.2025.